



# The most common and valuable accessories

Solution Simplified

## Efficient Functionality

Innovative Design

Outstanding Outcomes

Top-Quality Accessories



www.aczel.com

The Ultimate Answer for Your Applications

# The Most Common and Valuable Accessories

## Data Printer CPR 02



- GLP/GMP Header Footer.
  - Date & Time Password Protected.
  - Selectable Baud rate.
- \* Universal to connect any other brand.

## Anti-Vibration Pad VRD-03



Anti-vibration pad for analytical balance reduces vibrations.

## IQ, OQ, PQ Document LBA-X 20



validation processes used to ensure that an analytical balance is properly installed operates correctly, and performs within specified parameters.

## Anti-Theft Lock LBA-X 16



Provides Security and Prevents unauthorized access or theft of the Analytical Balance.

## Density Kit - CDK01



A density kit is a specialized accessory used with analytical balances to determine the density of solid and floating samples.

## Battery Backup LBA-X 01



A Battery backup ensures reliable power for an an during a power failure, delivering accurate results.

## 21 CFR PART 11 PC Software




21 CFR Part 11 is the part of Title 21 of the Code of Federal Regulations that establishes the United States Food and Drug Administration (FDA) regulations on electronic records and electronic signatures(ERES

## Specific Gravity Stand LBA-X 13.B



Below balance weighing is ideal for oversized or irregularly shaped objects, as well as hot/cold or magnetic materials that could interfere with the scale's accuracy.

## Test Weight



Calibration test weight ensures accurate balance calibration.

## Remote Display SRD - 01



Remote display for analytical balance provides convenient viewing from a distance

## Interface Cable LBA-X 15



Enables communication and data transfer between the balance and a computer or Data Printer

## Protection Cover LBA-X 03



Encloses the entire housing protecting the balances from stains and scratches