

Design

- Designed for easy operation, HS Series digital scales meet the needs of medical professionals and health practitioners worldwide.
- Balance standard of OIML III grade.



Application

- Hospital / Nursing Homes / Clinics
- Physicians
- Dietitian
- Telemedicine
- Gymnasiums / Health Clubs
- Shopping Malls
- Corporate Offices
- Defense
- Railway Stations / Airport
- Specially used in infant/baby weighing
- Ideal for hospital / medical / pediatric purpose
- Govt. Health Dept.
- Anganwadi (ICDS)

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

Construction

- Rugged & reliable design.
- Compact and space saving.

Features

- Large LCD display.
- Fast weighing response time.
- Auto power OFF.
- Auto height measured on display.
- Real time clock.

• Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)

BMI Height And Weight Scale

HS - BMI Series

Model	HS 200 BMI
Capacity (kg / lbs)	200 / 440
Readability (g)	100
Pan Size L x w (mm)	320 x 420 Anti-slip Platform Ensures Safety
Weighing Units	kg / lb
Max weighing (kg)	200
Division (g)	100
Stabilization Time (sec)	3
Battery	2 x 1.5V Battery
Height Rod (mm)	210 (Can Be Customized)
Height range to be measured	190cm (Optional 210cm)
Min value of height per division (mm)	5
Calibration	External
Display	LED Backlite / LCD
Power Supply	AC : 110/220V(±10%)
Operating Temperature	5° to 40°C
Net Weight (kg)	14
Gross Weight (kg)	15
Packing Size (L x W x H) (mm)	780 x 375 x 210