

Baby Scale

BS - IGMS Series



aczel®

...Since 1986



IGMS Series Pediatric Scale is designed as an all-in-one unit with no moving parts, making it a breeze to clean. The display results of a child comprising of all growth charts, as per World Health Organization (WHO). this system displays the CATEGORY of Underweight, Stunting and Wasting of a child.

Features

- USB port to copy the data on Pen Drive
- Data storage copy of min 500 Weighing.
- In Built Real Time Clock.
- Facility to entry ID, gender, age, height etc.
- Equipped with on fact / child growth monitoring system

Application :

- Specially used in infant/baby weighing
- Ideal for hospital / medical / pediatric purpose
- Govt. Health Dept.
- Anganwadi (ICDS)

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

• **Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)**

Model	BS 20 IGMS	BS 40 IGMS
Capacity (kg/lbs)	20 / 44	40 / 90
Accuracy (gm)	10	10
Function Keys	Music on Hold	Music on Hold
Pan Size (LxWxH)(mm)	550 x 270 x 240 with bowl	550 x 270 x 240 with bowl
Display	LCD Backlite	LCD Backlite
Power Adapter	6Vdc / 1amp	6Vdc / 1amp
Operation Temp.	10°C to 40°C (14°F to 104°F)	10°C to 40°C (14°F to 104°F)
Net Weight (kg)	4	4
Gross Weight (kg)	5	5
Packing Dimension (mm)	750 x 595 x 395	750 x 595 x 395