



The BS 20 Pediatric scale is designed as an all-in-one unit with no moving parts, making it a breeze to clean. Aczel averaging weight technology guarantees incredible sensitivity for your smallest patients. The BS 20 is lightweight, yet durable enough for the health care professional on the move. An easy to read bakelite LCD/bright LED makes this infant scale perfect for quick and accurate measurements. It also includes a self-adhesive measuring ruler (in fants - meter) to assess baby's height at the same time. Save Feature allows you to see what the last measurement was. Freeze features can hold the weight. All models feature Aczel topquality manufacturing for long life and owner satisfaction.

Features

- "AAA" Light weight battery operated scale.
- Convenient design for infant/baby weighing
- Tare and auto hold function
- Automatic zero setting
- Attachable baby length measure.
- Save & hold feature
- Display bright backlite LCD for easy viewing
- Response time max 3-5 sec.
- Horizontal leveling with height adjustable feet

Application :

- Specially used in infant/baby weighing
- Ideal for hospital / medical / pediatric purpose
- Govt. Health Dept.
- Anganwadi (ICDS)

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

• Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)

Baby Scale

BS Series

Model	BS 20
Capacity (kg/lbs)	20 / 44
Accuracy (gm)	5
Function Keys	Music on Hold
Over all Size (LxWxH)(mm)	545 x 330 x 145, ±5mm
Display	LCD Backlite
Power Adapter	6Vdc / 1amp
Operation Temp.	10°C to 40°C (14°F to 104°F)
Net Weight (kg)	4
Gross Weight (kg)	5
Packing Dimension (mm)	600 x 350 x 150