




Operation



Features

- High precision strain gauge sensors system.
- SS Cover or MS Powder Coated Body.
- Low battery and over load indication.
- Auto ON and OFF function.
- Step on technology calibration.
- Weight for Mom & Baby.
- Low Battery Indication

Application :

- Specially used in infant/baby weighing with Mother
- Ideal for hospital / medical / pediatric purpose
- Govt. Health Dept.
- Anganwadi (ICDS)

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- ISO 9001: 2015
- ISO 13485 : 2016

• Conforming to :- IS: 9281 (Part 1&2) / 1979

IS: 9281 (Part - 3) / 1981

IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)

Mother and Child Scale

HS Series

Model	HS 180
Capacity (kg / lbs)	180 / 400
Graduation (g)	d=100g
Display	Large LCD with backlite
Power	Lithium Battery(included)
Product Dimension (L x W x H)(mm)	300 x 300 x 28 ± 5 mm tol.
Net Weight (kgs) (±100 gm)	2
Gross Weight (kgs) (±100 gm)	2.2
Type of Material	Top Surface Stainless Steel

ACZET PVT. LTD.

Unit E2, Plot No. 15, WICEL Estate, Opp Seepz Gate No.1, Andheri (E), Mumbai 400 093.
Tel.: 022-4243 7748 • Mobile no: 8452002414 • Email: tender@aczel.com