



Specification :

- High Precision strain gauge sensor.
- LCD Display Size 3.2 inch.
- Low Battery / Over Load Indication.
- 6mm tempered safety glass platform with non slip base
- Energy saving automatic shutdown
- 3 Units Conversion : kg, lb, st.
- Automatic ON / OFF, Low Battery (LU) and Overload indication (Err)

Application :

- Specially used in infant/baby weighing
- Ideal for hospital / medical / pediatric purpose
- Govt. Health Dept.
- Anganwadi (ICDS)
- Health Club / GYM

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

• Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)

Health Scale

HS Series

Model	HS 180
Capacity (kg / lbs)	180 / 400
Accuracy (gm)	100
Pan Size (mm)	280 x 280 x 23.5
Display	Backlit : LCD / Size : 3.2 inch with back light & temperature
Power	2 x 1.5V (2AAA battery) / Lithium Battery
Net Weight (kg)	2
Gross Weight	2.5
Packing Dimension (L x W x H) (mm)	300 x 300 x 24