



### Design

- Designed for easy operation, HS Series digital scales meet the needs of medical professionals and health practitioners worldwide.
- Balance standard of OIML III grade.

### Construction

- Rugged & reliable design.
- Compact and space saving.

### Display

- Easy to read LED.

### Power Supply

- Power AC230V,  $\pm 10\%$ , 50-60 Hz.

### Features

- Built in 2 x 1.5V battery .
- Fast weighing response time.
- Auto power off.

## Height and Weight Scale

HS Series

Model	HS 200
Capacity	200 kg
Readability	100 g
Repeatability (S.D.)	0.1 kg
Linearity (+/-)	0.1 kg
Pan Size L x w (mm)	320 x 420 Anti-slip Platform Ensures Safety
Weighing Units	kg / LB
Max weighing	150/200kg
Division	100g
Stabilization Time (sec)	3
Battery	2 x 1.5V Battery
Height Rod (mm)	210 (Can Be Customized)
Height range to be measured	190cm (Optional 210cm)
Min value of height per division	5mm
Calibration	External
Display	LED Backlite LCD
Power Supply	AC : 110/220V( $\pm 10\%$ )
Operating Temperature	5° to 40°C
Net Weight / Gross Weight	14 / 15kg
CTN Meas. (mm)	780 x 375 x 210