



- Rugged & reliable design
- Compact and space saving.
- Easy to read VFD / backlit LCD display.
- Power AC 230V, + 10%, 50-60 Hz
- Built-in rechargeable battery.
- Piece counting function available.
- Optional RS 232C Interface.
- Auto power off.
- Auto Calibration.
- With windshield.
- Remote display.

Precision Balance

CG-M Series

Model	CG 602M
Capacity	600 gm
Readability	0.01 gm
Repeatability (+/-)	0.01 gm
Linearity (+/-)	0.02 gm
Pan Size (mm)	128 Ø