

Gold Testing Machine

STARK SDD PRO



- Detect 22 Metals at a time.
- Measuring results in 20 sec.
- Rapid, Accurate.

Model	Stark
Measuring Direction	Bottom to Top
Applied Application	Gold, Silver & Platinum Analysis in jewellery and alloys
X-ray Tube	Powerful 50 Watts air cooling X-ray tube
High Voltage	50kV (1.2mA) Software Control Optimized
Detector	High performance SDD
Measurement Time	10 - 40 sec.
Collimator	0.3mmØ or 0.5mmØ, Optional four positions Collimator changer
Sample Stage	Fixed Sample Positioning
Power Supply	220ACV 50/60HZ
Inside Chamber	300 x 310 x 104 mm (WxDxH)
Dimension	400 x 405 x 455 mm (WxDxH)
Weight	40 kg

Specifications are subject to change without prior notice.
Images shown are for illustration purpose only.