



Standard Features

- High Resolution up to 30,000Conts.
- Built-in rechargeable battery for continuous use.
- Piece counting function available.
- Environmentally protected Load cell.
- Rust Proof Stainless Steel Platter.
- Auto Calibration.
- Stainless Steel Platter & Body.
- RS232C
- Piece counting facility

Application

- Engineering Industry
- Ginning industry
- Grain Industry
- Laboratory Industry
- Chemical Industry
- Pharma Industry & others

Table Top Scale

CG -S Series

Model	CG 3S	CG 6S	CG 10S	CG 15S	CG 30S
Maximum Capacity	3 kg	6 kg	10 kg	15 kg	30 kg
Readability	0.1 gm	0.5 gm	0.5 gm	1 gm	1 gm
Pan Size (mm)	260 x 210 mm				
Response Time	2 - 3 Sec.				
Display	LCD Display				
Calibration	External Calibration				
Operating Temp.	5° to 40°C				
Power Supply	AC 230V, +/- 10%, 50-60Hz				

Typical minimum sample weight according to USP Chapter41