



### Specification

- LCD display size: 102.5mm\*67.5mm with backlight
- Auto switch
- Auto zero, zero off
- Low battery /overload indication
- You can choose any colors
- Body weight/ Fat/ Hydration/ Muscle/Bones/ Calories/ BMI with ultrasonic height function

### Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

- Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)

## Ultrasonic Height Body Fat Analyzer

HS - PH Series

| Model                             | HS 200PH   |
|-----------------------------------|--|
| Capacity (kg / lbs)               | 200 / 440  |
| Accuracy (g)                      | 50   |
| Display Type                      | LCD  |
| Unit                              | Kg, LB, ST   |
| Function                          | Body Fat / Water / Muscle / Bones / Calories / BMI / Height / Weight |
| Power                             | 4 x 1.5v 7 AAA battery (included) or the 6~12v dc power supply       |
| Scale Dimensions (L x W x H) (mm) | 448 x 340 x 2075   |
| Net weight (kg)                   | 4.55   |
| Gross weight (kg)                 | 5.93   |