



**CUB 5**



**CUB 10 ~ 30**

ACZET Ultrasonic units are used in Industrial and Analytical Laboratories. With their high power, optimized parameter settings and strict parameter monitoring and assisted by a large range of useful accessories the ultrasonic units guarantee safe and reliable results for the multifaceted laboratory applications. The ACZET Laboratory equipment is complete with a full programme of cleaning chemicals developed and produced in our own laboratory.

**Application :  
Laboratory**

Get consistent cleaning with every cycle. Our Ultrasonic baths thoroughly remove blood, protein, and contaminants. Use it for everything from glassware to lenses to instruments and precision components. Beyond Cleaning ACZET ultrasonic cleaners are extremely versatile. They can also be used in laboratory applications for cell separation, mixing emulsifying, sample preparation and degassing of liquids.

**Electronics**

Get a safe, unwavering cleaning for even the most delicate of components. Aczet baths completely remove flux and contaminants from such precision parts as PC boards, SMDs, quartz crystals, capacitors, and many others.

**Jewellery**

See the difference and performance of Aczet ultrasonic baths. Thoroughly clean and restore brilliance to watches, chains and charms, settings, coins, fine Industrial The demands of light industrial workflow require an ultrasonic bath that can stand up to the challenge. Our baths are perfect for deep cleaning to remove dirt, grease, waxes, and oils from light industrial parts and components of all kinds - including steel, light & nonferrous metals, plastic and glass. jewellery and clockworks of every type.

Model	CUB - 2.5	CUB - 5	CUB - 10	CUB - 20	CUB - 30
Capacity	2.5 Ltr.	4.5 Ltr.	9.8 Ltr.	20 Ltr.	27 Ltr.
Input	230 Vac, ±20%, 50/60Hz				
Frequency	40KHz Ultrasonic power 120W				
Adjustable Time	1 to 99 minutes	1 to 99 minutes	1 to 99 minutes	1 to 99 minutes	1 to 99 minutes
Adjustable Temperature	0°C to 80°C	0°C to 80°C	0°C to 80°C	20°C to 80°C	20°C to 80°C
Heater Power	50 W	120W	240W	600W	600W
Tank Size (mm)	150 x 137 x 150	300 x 155 x 100	300 x 240 x 150	500 x 300 x 150	500 x 300 x 200
Packing Size (mm)	250 x 200 x 250	350 x 250 x 330	440 x 405 x 390	640 x 450 x 390	640 x 450 x 390
Gross Weight	4.2 kg	6 kg	11 kg	21 kg	31 kg