



# Animal Weighing Scale

CTG - A Series



The Veterinary Scale offers a convenient, compact and portable weighing capability with high capacity. This durable product offers unparalleled weighing performance, with the advantages of low profile and affordability. Standard features include lbs, kg and weighing modes, weight hold feature, stable weight indication and auto power off to extend battery life. It's Ideal for use in veterinary clinics, kennels, or grooming shops. Removable stainless steel platform cleans up easily. Battery powered lighted display retains weight after animal is removed from scale.

Our exclusive MotionLock feature allows an accurate weight to be held regardless of user motion. The scale software was specifically designed to provide an accurate weight even if the patient is unsteady and is unable to remain perfectly still.



**Features**

- Weighing units (kg, lb, pcs)
- Aluminum case indicator
- Hold function
- Full range tare
- Zero Tracking
- Tare function
- Wall mount bracket
- Large stainless steel platform
- Modern low profile design
- Simple 4 button operation
- USB or RS-232 interface
- Auto sleep / power down function to save battery life
- Low battery indication
- AC Adapter

**Approval :**

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

• **Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)**

Model	CTG 150A	CTG 200A
Capacity (kg / lbs)	150 / 330	200 / 440
Readability (g)	50	50
Display	Large backlit LCD display	
Dimensions (mm)	600 x 500	900 x 600