

Handrail Scale

HS Series



Application

- Hospital / Nursing Homes / Clinics
- Physicians
- Dietitian
- Telemedicine
- Gymnasiums / Health Clubs
- Shopping Malls
- Corporate Offices
- Defense
- Railway Stations / Airport

This Handrail Scale is ideal for extra patient support during weighing. The large durability. The maximum capacity is up to 300kg

- **Power Supply : AC Adaptor & Rechargeable Battery**

Approval :

- Models Approved by Legal Metrology (WM). Govt. of India.
- BIS Certified ISI - 7800029512
- ISO 9001: 2015
- ISO 13485 : 2016

- **Conforming to :- IS: 9281 (Part 1&2) / 1979, IS: 9281 (Part - 3) / 1981, IS: 9281 (Part - 4) / 1983 (Reaffirmed 2006)**

Model	HS 300
Capacity (kg / lbs)	300 / 660
Accuracy (g)	100
Net Weight Approx. (kg)	16.5
Electronic Supply	230V/110V, 50Hz
Platform Size (mm)	400 x 450 x 95